

May 20, 2024 (Ver.2)
Approval No. 22900BZX00291000
Class IV Artificial Dura Mater(70511000)
DuraBeam®
 Sterilized, Single Use
 Thickness : 0.30mm



Sizes: LL (120mmx120mm, D-300-1212)
 L (100mmx100mm, D-300-1010)
 M (100mmx50mm, D-300-1005)
 S (50mmx50mm, D-300-0505)

"DuraBeam®" Synthetic Artificial Dura Mater Ion-irradiated ePTFE sheet Usage Warning

- **DuraBeam®** synthetic artificial dura mater has a Front Surface (processed by ion beam irradiation) and a Back Surface (unprocessed).
 (The material is expanded polytetrafluoroethylene (ePTFE), which is highly biocompatible, safe and stable.)

Fig.1 DuraBeam® Front Surface (Ion-beam irradiated)

Fig.2 DuraBeam® Back Surface (unprocessed)



- "ION" or "Ion" engraved
- Matte white to very slightly brown

- The engraving is reverse
- Smooth and glossy white

- Please use **DuraBeam®** with the Front Surface (the ion-irradiated surface) facing the biological dura mater. Do not apply **DuraBeam®** directly to brain tissue. For details, please refer to the attached document.

Fig.3 Standard/Recommended Usage
 Implant inlay under dura mater, with the Front Surface (the ion-irradiated surface) facing outwards.





Fig.4 Avoid unless there is a special reason
 We do not recommend "onlay", placing it over the dura mater, and implanting it with the Front Surface (the ion-irradiated surface) facing inwards. (When performing the Figure 4 usage, please be careful not to let DuraBeam® come into contact with brain tissue.)



It has been reported that the biological dura mater adheres and regenerates onto the Front Surface of DuraBeam®, reducing the adverse event of cerebrospinal fluid leakage.